

## Specialty Items

### *Finneran Products Certified For Science™*

#### Capsule Holders



For use in dissolution vessels. Three prongs hold the capsule in place, permitting the completion of dissolution tests without capsule flotation.

**Sold individually**

Cat. No.	Description
91700-00	Black Capsule Holder
91700-01	Black Metallic Capsule Holder

#### Ampoule Holders



The Finneran Ampoule Storage Holder seals itself to your open ampule, allowing for easier sample storage, and eliminating the need to transfer the contents.

The added spring allows for easy access and removal of the ampule.

Just break off the 2mL ampule top, place the ampule in our storage holder, and attach the closure. Easily inject through the open top screw cap/stopper to remove your sample.

Made from K-Resin

**Sold individually**

Cat. No.	Description
9450-01	2mL Storage Ampoule Holder with Spring, Stopper and White Polypropylene Closure for Snapped Open Ampoules